

Hassle Free Australia Wide Delivery



Assembled and Ready to Drive

All you have to do is Lift the seat back, stand up the tiller and swap the arm rests from side to side.

Turn the key and drive you new scooter off the platform. Your new scooter has been through a Pre Delivery Service and Testing. The Batteries are fully charged.



Slide the seat back. The slider flap is just under the right edge of the seat. Lift up the tiller. Push in the locator lever at the bottom of the tiller. Loosen the hand screw at the back of the seat and simply swap the armrests side to side and make upright. Its that easy. You are now Mobile

Easy to Un Pack

Its like an upside down shoe box. Use a Sharp Knife to remove the plastic wrapping. Cut the blue straps and lift the box off your new scooter. Remove the foam and cardboard packaging.



Securely Packaged and Wrapped

First class packaging to ensure your new scooter arrives exactly as it left the factory Plastic wrapped to survive any circumstance



Tailgate loader to you door

or Pick up from the local depot

Your new scooter is shipped on a skip pallet for easy handling and delivery. The freight company can use a tailgate loader to deliver the scooter to your garage or yard. In some places it is easier and cheaper to simply go to the depot and drive your scooter home.



4,000 plus deliveries and No complaints so far